

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 012 BL USP

Vol. flask, 10 ml, NS 10, blue grad., USP, A, conf. cert.

Batch: 31-21


Serial:

Measuring Nr. 48405

Batch size:

500

| | |
|------------------------|-----------------|
| Volume: | 10,0000 ml |
| Tolerance: | $\pm 0,0200$ ml |
| Average value: | 10,0031 ml |
| Deviation: | 0,0043 ml |
| Precision (CV): | 0,0427 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

17.08.2021

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M