

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 020 BL

**Vol. flask, 200 ml, NS 14, blue grad., A, conf. cert.**


**Batch:** 51-20

**Serial:**

**Measuring Nr.** 43547

**Batch size:** 1198

|                        |                 |
|------------------------|-----------------|
| <b>Volume:</b>         | 200,000 ml      |
| <b>Tolerance:</b>      | $\pm 0,1500$ ml |
| <b>Average value:</b>  | 200,0022 ml     |
| <b>Deviation:</b>      | 0,0243 ml       |
| <b>Precision (CV):</b> | 0,0121 %        |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                          |
|------------------|--------------------------|
| <b>Operator:</b> | 7                        |
| <b>Date:</b>     | 13.01.2021               |
| <b>Balance:</b>  | Sartorius - LA<br>1200 S |

DE-M