

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 025 BL USP

Vol. flasks, 250 ml, NS 14, blue, USP, A, conf. cert.


Batch: 09-22

Serial:

Measuring Nr. 53111

Batch size: 276

| | |
|------------------------|-----------------|
| Volume: | 250,000 ml |
| Tolerance: | $\pm 0,1200$ ml |
| Average value: | 250,0083 ml |
| Deviation: | 0,0213 ml |
| Precision (CV): | 0,0085 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 17.03.2022

Balance: Sartorius - LA
1200 S

DE-M