

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 025 BL

**Vol. flask, 250 ml, NS 14, blue grad., A, conf. cert.**


**Batch:** 26-20

**Serial:**

**Measuring Nr.** 39132

**Batch size:** 1177

|                        |                   |
|------------------------|-------------------|
| <b>Volume:</b>         | <b>250,000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,1500$ ml   |
| <b>Average value:</b>  | 250,0073 ml       |
| <b>Deviation:</b>      | 0,0169 ml         |
| <b>Precision (CV):</b> | 0,0068 %          |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                          |
|------------------|--------------------------|
| <b>Operator:</b> | 7                        |
| <b>Date:</b>     | 21.07.2020               |
| <b>Balance:</b>  | Sartorius - LA<br>1200 S |

DE-M