

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 662 025 BL

Vol. flask, 250 ml, NS 14, blue grad., A, conf. cert.

Batch: 25-24


Serial:

Measuring Nr. 67846

Batch size:

496

Volume:	250,000 ml
Tolerance:	$\pm 0,1500$ ml
Average value:	250,0034 ml
Deviation:	0,0194 ml
Precision (CV):	0,0078 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:

7

Date:

04.07.2024

Balance:

Sartorius - LA  
1200 S

DE-M