


# Certificate of performance

Article number 3 662 025 USP  
 Batch 02-17  
 Serial number  
 Measuringnr./Batch size 10054 / 200  
 Balance Sartorius - LA 1200 S  
 Operator 4  
 Date 19.01.2017

## Volumetric Flask 250 ml Class A NS 14 USP PE-plug

Volume 250,0000ml  
 Tolerance 0,1200ml  
 Average volume 249,9836 ml  
 Deviation 0,0198 ml  
 Number of measurings 15



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
 97877 Wertheim  
 Germany  
 +49 (0) 9342 / 9301-0  
 | info@witeg.de

