

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 040

Volumetric flask, 400 ml, NS 19, A, KB

Batch: 14-24


Serial:

Measuring Nr. 66784

Batch size:

90

| | |
|------------------------|-------------------|
| Volume: | 400,000 ml |
| Tolerance: | $\pm 0,2500$ ml |
| Average value: | 399,9672 ml |
| Deviation: | 0,0387 ml |
| Precision (CV): | 0,0097 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|--------------------------|
| Operator: | 4 |
| Date: | 16.04.2024 |
| Balance: | Sartorius - LA 1200 S |

DE-M