

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 060 BL

Volumetric flask, 600 ml, blue, NS 24, A, conf.cert.

Batch: 36-24


Serial:

Measuring Nr. 68941

Batch size:

73

Volume:	600,000 ml
Tolerance:	$\pm 0,4000$ ml
Average value:	600,1118 ml
Deviation:	0,0279 ml
Precision (CV):	0,0047 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

05.09.2024

Balance:

Sartorius - LA
1200 S

DE-M