

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 101 BL

Vol. flasks, 1000 ml, NS 29, blue grad., A, conf. cert.


Batch: 16-18

Serial:

Measuring Nr. 20362

Batch size: 65

| | |
|------------------------|---------------------|
| Volume: | 1000,0000 ml |
| Tolerance: | $\pm 0,6000$ ml |
| Average value: | 999,9553 ml |
| Deviation: | 0,0600 ml |
| Precision (CV): | 0,0060 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 44

Date: 07.05.2018

Balance: Sartorius -
LE244S-OCE

DE-M