

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 200 BL USP

Vol. flask, 2000 ml, NS 29, blue grad., USP, A, conf. cert


Batch: 45-18

Serial:

Measuring Nr. 25657

Batch size: 50

Volume:	2000,0000 ml
Tolerance:	± 0,5000 ml
Average value:	2000,1276 ml
Deviation:	0,0417 ml
Precision (CV):	0,0021 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	17.11.2018
Balance:	Sartorius - LC 6200 S

DE-M