

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 662 300

Volumetric flasks, 3000 ml, NS 29, A, conf. cert.


Batch: 30-18

Serial:

Measuring Nr. 23342

Batch size: 20

Volume:	3000,0000 ml
Tolerance:	$\pm 0,8000$ ml
Average value:	3000,1834 ml
Deviation:	0,0460 ml
Precision (CV):	0,0015 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	14.08.2018
Balance:	Sartorius - LE244S-OCE

DE-M