

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 300

Volumetric flask, 3000 ml, NS 29, A, conf. cert.


Batch: 07-23

Serial:

Measuring Nr. 61807

Batch size: 25

Volume:	3000,0000 ml
Tolerance:	± 0,8000 ml
Average value:	3000,1555 ml
Deviation:	0,1285 ml
Precision (CV):	0,0043 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	23.02.2023
Balance:	Sartorius - MSE14202S- OCE-D0

DE-M