

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 400 BL

**Vol. flasks, 4000 ml, NS 34, blue grad., A, conf. cert.**

**Batch:** 02-18


**Serial:**

**Measuring Nr.** 16932

**Batch size:**

6

<b>Volume:</b>	4000,0000 ml
<b>Tolerance:</b>	$\pm 1,2000$ ml
<b>Average value:</b>	3999,8666 ml
<b>Deviation:</b>	0,1516 ml
<b>Precision (CV):</b>	0,0038 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	15.01.2018
<b>Balance:</b>	Sartorius - LC 6200 S

DE-M