

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 662 400 BL

Vol. flask, 4000 ml, NS 34, blue grad., A, conf. cert.


Batch: 14-24

Serial:

Measuring Nr. 67070

Batch size: 15

Volume:	4000,0000 ml
Tolerance:	$\pm 1,2000$ ml
Average value:	4000,5127 ml
Deviation:	0,2397 ml
Precision (CV):	0,0060 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	03.05.2024
Balance:	Sartorius - MSE14202S- OCE-DO

DE-M