

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 600

Volumetric flask, 6000 ml, NS 34, A, conf.cert.


Batch: 14-24

Serial:

Measuring Nr. 66540

Batch size: 13

Volume:	6000,000 ml
Tolerance:	$\pm 1,2000$ ml
Average value:	6000,5000 ml
Deviation:	0,1459 ml
Precision (CV):	0,0024 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	03.04.2024
Balance:	Sartorius - MSE14202S- OCE-DO

DE-M