

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 000 5

Vol. flasks, 0,5 ml, NS 7, trapezoidal, A, conf. cert.

Batch: 30-17


Serial:

Measuring Nr. 12640

Batch size:

50

| | |
|------------------------|------------------|
| Volume: | 0,5000 ml |
| Tolerance: | $\pm 0,0250$ ml |
| Average value: | 0,5021 ml |
| Deviation: | 0,0023 ml |
| Precision (CV): | 0,4588 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 9 |
| Date: | 03.08.2017 |
| Balance: | Sartorius - LE244S-OCE |

DE-M