

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 000

Vol. flasks, 1 ml, NS 7, trapezoidal, A, conf. cert.

Batch: 01-18


Serial:

Measuring Nr. 16854

Batch size:

500

Volume:	1,0000 ml
Tolerance:	± 0,0250 ml
Average value:	1,0013 ml
Deviation:	0,0030 ml
Precision (CV):	0,2964 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 9

Date: 09.01.2018

Balance: Sartorius -
LE244S-OCE

DE-M