

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 00010

Volumetric flask, 1 ml, NS 10, trapezoidal, A, conf. cert.


Batch: 49-19

Serial:

Measuring Nr. 34893

Batch size: 90

Volume:	1,0000 ml
Tolerance:	± 0,0400 ml
Average value:	0,9930 ml
Deviation:	0,0063 ml
Precision (CV):	0,6277 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	04.12.2019
Balance:	Sartorius - LE244S-0CE

DE-M