

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 001

Vol. flasks, 2 ml, NS 7, trapezoidal, A, conf. cert.


Batch: 28-18

Serial:

Measuring Nr. 22107

Batch size: 1100

Volume:	2,0000 ml
Tolerance:	± 0,0250 ml
Average value:	2,0017 ml
Deviation:	0,0050 ml
Precision (CV):	0,2476 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 44

Date: 10.07.2018

Balance: Sartorius -
LE244S-OCE

DE-M