

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 015 USP

Volumetric flask, 5 ml, USP, NS 10, trapez., A, conf. cer


Batch: 18-20

Serial:

Measuring Nr. 37526

Batch size: 112

Volume:	5,0000 ml
Tolerance:	$\pm 0,0800$ ml
Average value:	5,0084 ml
Deviation:	0,0041 ml
Precision (CV):	0,0816 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	28.04.2020
Balance:	Sartorius - LE244S-OCE

DE-M