

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 050 BL

Vol. flask, 50 ml, NS 12, trapez., blue, A, conf. cert.


Batch: 06-19

Serial:

Measuring Nr. 27958

Batch size: 100

Volume:	50,0000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	49,9840 ml
Deviation:	0,0101 ml
Precision (CV):	0,0202 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 4

Date: 25.02.2019

Balance: Sartorius -
LE244S-OCE

DE-M