

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 002

Vol. flask, 20 ml, NS 10, amber glass, A, conf. cert.

Batch: 32-19


Serial:

Measuring Nr. 32782

Batch size:

1200

Volume:	20,0000 ml
Tolerance:	± 0,0400 ml
Average value:	20,0045 ml
Deviation:	0,0090 ml
Precision (CV):	0,0448 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

05.09.2019

Balance:

Sartorius -
LE244S-OCE

DE-M