

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 005

Vol. flasks, 50 ml, NS 12, amber glass, A, conf. cert.

Batch: 30-18


Serial:

Measuring Nr. 23111

Batch size:

2620

Volume:	50,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	50,0099 ml
Deviation:	0,0140 ml
Precision (CV):	0,0280 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	08.08.2018
Balance:	Sartorius - LE244S-OCE

DE-M