


# Certificate of performance

Article number 3 670 011  
 Batch 05-17  
 Serial number  
 Measuringnr./Batch size 10081 / 500  
 Balance Sartorius - LE244S-0CE  
 Operator 4  
 Date 08.02.2017

## Volumetric flask 100 ml Ciss A NS 14 KB brown glass

Volume 100,0000ml  
 Tolerance 0,1000ml  
 Average volume 100,0096 ml  
 Deviation 0,0141 ml  
 Number of measurings 20



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
 97877 Wertheim  
 Germany  
 +49 (0) 9342 / 9301-0  
 | info@witeg.de

