

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 012 USP

Vol. flasks, 10 ml, NS 10, amber glass, A, conf. cert.


Batch: 46-21

Serial:

Measuring Nr. 51177

Batch size: 499

| | |
|------------------------|--------------------|
| Volume: | 10,0000 ml |
| Tolerance: | ± 0,0200 ml |
| Average value: | 9,9968 ml |
| Deviation: | 0,0071 ml |
| Precision (CV): | 0,0715 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 17.12.2021 |
| Balance: | Sartorius - LE244S-OCE |

DE-M