


Certificate of performance

Article number 3 670 012
Batch 05-17
Serial number
Measuringnr./Batch size 10079 / 500
Balance Sartorius - LE244S-0CE
Operator 4
Date 06.02.2017

Volumetric flask 10 ml Class A NS 10 brown glass

Volume 10,0000ml
Tolerance 0,0400ml
Average volume 9,9943 ml
Deviation 0,0069 ml
Number of measurings 20



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
| info@witeg.de



ISO 9001
CERTIFIED

