

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 670 013 USP

**Vol. flask, 20 ml, USP, NS 12, amber gl., A, conf. cert.**

**Batch:** 26-19


**Serial:**

**Measuring Nr.** 31157

**Batch size:**

141

<b>Volume:</b>	20,000 ml
<b>Tolerance:</b>	$\pm 0,0800$ ml
<b>Average value:</b>	20,0012 ml
<b>Deviation:</b>	0,0077 ml
<b>Precision (CV):</b>	0,0386 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

01.07.2019

**Balance:**

Sartorius -  
LE244S-OCE

DE-M