

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 013

Vol. flasks, 20 ml, NS 12, amber glass, A, conf. cert.


Batch: 04-21

Serial:

Measuring Nr. 44253

Batch size: 498

| | |
|------------------------|-----------------|
| Volume: | 20,000 ml |
| Tolerance: | $\pm 0,0600$ ml |
| Average value: | 20,0016 ml |
| Deviation: | 0,0082 ml |
| Precision (CV): | 0,0410 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 05.02.2021 |
| Balance: | Sartorius - LE244S-OCE |

DE-M