

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 017 PE

Volumetric


Batch: 12-19

Serial:

Measuring Nr. 28640

Batch size: 298

Volume:	20,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	19,9929 ml
Deviation:	0,0089 ml
Precision (CV):	0,0443 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	0
Date:	23.03.2019
Balance:	Sartorius - LE244S-OCE

