

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 670 017

Vol. flasks, 20 ml, NS 14, amber glass, A, conf. cert.


Batch: 50-22

Serial:

Measuring Nr. 60728

Batch size: 280

Volume:	20,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	19,9885 ml
Deviation:	0,0136 ml
Precision (CV):	0,0680 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	14.12.2022
Balance:	Sartorius - MSE224S-0CE -DU

