

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 670 018 P

**Vol. flasks, 25 ml, NS 14, Amber glass, A, conf. cert.**

**Batch:** 07-18

**Serial:**


-

**Measuring Nr.** 18991

**Batch size:**

280

<b>Volume:</b>	<b>25,0000 ml</b>
<b>Tolerance:</b>	$\pm 0,0600$ ml
<b>Average value:</b>	24,9595 ml
<b>Deviation:</b>	0,0060 ml
<b>Precision (CV):</b>	0,0240 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

44

**Date:**

19.03.2018

**Balance:**

Sartorius -  
LE244S-OCE

DE-M