

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 670 018

Vol. flasks, 25 ml, NS 14, amber glass, A, conf. cert.

Batch: 07-18


Serial: -

Measuring Nr. 18991

Batch size:

280

Volume:	25,0000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	24,9595 ml
Deviation:	0,0060 ml
Precision (CV):	0,0240 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:

44

Date:

19.03.2018

Balance:

Sartorius -  
LE244S-0CE

DE-M