

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 020 USP

Vol. flask, 200 ml, NS 14, amber glass, A, conf. cert.

Batch: 35-21


Serial:

Measuring Nr. 48809

Batch size:

278

Volume:	200,000 ml
Tolerance:	$\pm 0,1000$ ml
Average value:	199,9691 ml
Deviation:	0,0239 ml
Precision (CV):	0,0120 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

13.09.2021

Balance:

Sartorius -
MSE14202S-
OCE-D0

DE-M