

Certificate of performance

Article: 3 670 025 19 USP

Vol. flask, 250 ml, NS 19, amber glass, A, conf. cert.


Batch: 02-19

Serial:

Measuring Nr. 26776

Batch size: 90

Volume:	250,000 ml
Tolerance:	$\pm 0,2000$ ml
Average value:	249,9904 ml
Deviation:	0,0272 ml
Precision (CV):	0,0109 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

14.01.2019

Balance:

Sartorius - LA
1200 S