

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 025 USP

Vol. flask, 250 ml, Ns 14, amber gl., USP, A, conf.cert.

Batch: 48-20


Serial:

Measuring Nr. 43087

Batch size:

148

Volume:	250,0000 ml
Tolerance:	$\pm 0,1200$ ml
Average value:	249,9854 ml
Deviation:	0,0262 ml
Precision (CV):	0,0105 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

08.12.2020

Balance:

Sartorius - LA
1200 S

DE-M