

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 101

Vol. flasks, 1000 ml, NS 29, amber glass, A, conf. cert.


Batch: 18-18

Serial:

Measuring Nr. 20792

Batch size: 61

Volume:	1000,000 ml
Tolerance:	$\pm 0,6000$ ml
Average value:	999,9961 ml
Deviation:	0,0742 ml
Precision (CV):	0,0074 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 4

Date: 17.05.2018

Balance: Sartorius - LC
6200 S

DE-M