

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 200 USP

Vol. flask, 2000 ml, NS 29, amber glass, A, conf. cert.


Batch: 43-18

Serial:

Measuring Nr. 25658

Batch size: 21

Volume:	2000,0000 ml
Tolerance:	$\pm 0,5000$ ml
Average value:	2000,0554 ml
Deviation:	0,0714 ml
Precision (CV):	0,0036 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 4

Date: 17.11.2018

Balance: Sartorius - LC
6200 S

DE-M