

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 674 0760

Control volumetric flask, 740-760 ml, A, conf. cert.




Batch: 18-21

Serial:

Measuring Nr. 46553

Batch size: 5

Volume:	740,000 ml	750,000 ml	760,000 ml
Tolerance:	± 0,5000 ml	± 0,4000 ml	± 0,5000 ml
Average value:	740,2317 ml	750,1933 ml	760,1506 ml
Deviation:	0,0702 ml	0,0938 ml	0,0929 ml
Precision (CV):	0,0095 %	0,0125 %	0,0122 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	07.05.2021
Balance:	Sartorius - LA 1200 S

DE-M