

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 420 080

Vol. pipette f. determin. of milk, 10,75 ml, conf. cert.


Batch: 25-19

Serial: -

Measuring Nr. 31148

Batch size: 5

Volume:	10,7500 ml
Tolerance:	$\pm 0,0200$ ml
Average value:	10,7542 ml
Deviation:	0,0018 ml
Precision (CV):	0,0163 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to -

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

28.06.2019

Balance:

Sartorius -
LE244S-OCE

DE-M