

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 011

Grad. pipettes, 0,1:0,01 ml, AS, conf. cert.




Batch: 05-18

Serial:

Measuring Nr. 18312

Batch size: 280

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,0984 ml	0,0583 ml	0,0178 ml
Deviation:	0,0015 ml	0,0009 ml	0,0009 ml
Precision (CV):	1,4710 %	1,5275 %	4,5991 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	21.02.2018
Balance:	Sartorius - LE244S-OCE

DE-M