

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 011

Grad. pipettes, 0,1:0,01 ml, AS, conf. cert.

Batch: 31-21




Serial:

Measuring Nr. 48140

Batch size:

276

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,0992 ml	0,0593 ml	0,0183 ml
Deviation:	0,0019 ml	0,0013 ml	0,0007 ml
Precision (CV):	1,8907 %	2,2320 %	3,4344 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

03.08.2021

Balance:

Sartorius -
LE244S-OCE

DE-M