

Certificate of performance



Article: 3 510 011

Grad. pipettes, 0,1:0,01 ml, AS, conf. cert.




Batch: 09-23

Serial:

Measuring Nr. 61878

Batch size: 280

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,1006 ml	0,0605 ml	0,0193 ml
Deviation:	0,0014 ml	0,0005 ml	0,0004 ml
Precision (CV):	1,3574 %	0,7936 %	2,1963 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	06.03.2023
Balance:	Sartorius - LE244S-OCE

