






Qualitätszertifikat / Certificate of performance

Cat. No. 3.510.021 BL	Graduated pipettes, measure AS	Messpipetten AS	
Batch No. 15-17		Code: 83333	
Testvolumen Test volume:	0.2000 ml	Variationskoeffizient Coefficient of variation:	0.6525 %
Toleranz Tolerance:	0.0060 ml	Richtigkeit Accuracy:	-0.4550 %
			
Testvolumen Test volume:	0.1200 ml	Variationskoeffizient Coefficient of variation:	1.2199 %
Toleranz Tolerance:	0.0060 ml	Richtigkeit Accuracy:	-1.5417 %
			
Testvolumen Test volume:	0.0400 ml	Variationskoeffizient Coefficient of variation:	2.0147 %
Toleranz Tolerance:	0.0060 ml	Richtigkeit Accuracy:	-5.6250 %
			
<small> Diese Angaben sind Endprüfwerte bezogen auf 20°C, d.h. gleiche Temperatur von Gerät und Umgebung und Flüssigkeit(aqua dest.). Die Prüfung erfolgte gemäß ISO 835 These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20°C). The test was effected according to ISO 835 </small>			
Waage / Balance: Sartorius LE244S-0CE		Datum / Date: 26.04.2017	
Prüfer / Operator: 9			
	Witeg Labortechnik GmbH Am Bildacker 16 97877 Wertheim Tel: +49 9342 93010 http://www.witeg.de		ISO 9001 CERTIFIED