

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 021

Grad. pipettes, 0,2:0,01 ml, AS, conf. cert.




Batch: 26-18

Serial:

Measuring Nr. 21867

Batch size: 280

Volume:	0,2000 ml	0,1200 ml	0,0400 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,1993 ml	0,1185 ml	0,0390 ml
Deviation:	0,0023 ml	0,0008 ml	0,0006 ml
Precision (CV):	1,1452 %	0,6811 %	1,4429 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 44

Date: 03.07.2018

Balance: Sartorius -
LE244S-OCE

DE-M