

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 021

**Grad. pipettes, 0,2:0,01 ml, AS, conf. cert.**




**Batch:** 07-20

**Serial:**

**Measuring Nr.** 36209

**Batch size:** 280

Volume:	0,2000 ml	0,1200 ml	0,0400 ml
<b>Tolerance:</b>	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
<b>Average value:</b>	0,1991 ml	0,1204 ml	0,0392 ml
<b>Deviation:</b>	0,0020 ml	0,0015 ml	0,0005 ml
<b>Precision (CV):</b>	0,9867 %	1,2905 %	1,3588 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	17.02.2020
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

