

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 021

Grad. pipettes, 0,2:0,01 ml, AS, conf. cert.

Batch: 44-21




Serial:

Measuring Nr. 50271

Batch size:

280

Volume:	0,2000 ml	0,1200 ml	0,0400 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,2004 ml	0,1214 ml	0,0398 ml
Deviation:	0,0016 ml	0,0011 ml	0,0009 ml
Precision (CV):	0,7976 %	0,9011 %	2,2192 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

08.11.2021

Balance:

Sartorius -
LE244S-OCE

DE-M