

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 101 BL

**Graduated pipettes, 1 ml, AS, conf. cert.**

**Batch:** 26-17




**Serial:**

**Measuring Nr.** 11974

**Batch size:**

500

| Volume:         | 1,0000 ml   | 0,6000 ml   | 0,2000 ml   |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value:  | 0,9993 ml   | 0,5972 ml   | 0,1983 ml   |
| Deviation:      | 0,0015 ml   | 0,0015 ml   | 0,0011 ml   |
| Precision (CV): | 0,1501 %    | 0,2430 %    | 0,5501 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                           |
|-----------|---------------------------|
| Operator: | 44                        |
| Date:     | 10.07.2017                |
| Balance:  | Sartorius -<br>LE244S-0CE |

