

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 101 BL

Grad. pipettes, 1:0,01 ml, blue grad., AS, conf. cert.




Batch: 02-24

Serial:

Measuring Nr. 65759

Batch size: 499

| Volume: | 1,0000 ml | 0,6000 ml | 0,2000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0070 ml | ± 0,0070 ml | ± 0,0070 ml |
| Average value: | 0,9977 ml | 0,6013 ml | 0,2007 ml |
| Deviation: | 0,0015 ml | 0,0013 ml | 0,0008 ml |
| Precision (CV): | 0,1487 % | 0,2226 % | 0,3878 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 16.01.2024 |
| Balance: | Sartorius - LE244S-OCE |

