

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 101 USP

**Grad. pipettes, 1:0,01 ml, AS, USP, conf. cert.**




**Batch:** 21-21

**Serial:**

**Measuring Nr.:** 47115

**Batch size:** 168

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
<b>Tolerance:</b>	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
<b>Average value:</b>	0,9995 ml	0,5972 ml	0,1994 ml
<b>Deviation:</b>	0,0019 ml	0,0012 ml	0,0013 ml
<b>Precision (CV):</b>	0,1867 %	0,2025 %	0,6317 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	09.06.2021
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

