

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 101 USP

**Grad. pipettes, 1:0,01 ml, AS, USP, conf. cert.**




**Batch:** 38-21

**Serial:**

**Measuring Nr.** 49416

**Batch size:** 227

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	1,0005 ml	0,6018 ml	0,2010 ml
Deviation:	0,0026 ml	0,0013 ml	0,0009 ml
Precision (CV):	0,2577 %	0,2239 %	0,4747 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	30.09.2021
Balance:	Sartorius - LE244S-OCE

